micros *Austria*

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
EN

Microtome blades

Disposable blades for sliding and rotary microtome

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www.micros.at

CH REP

Biosystems Switzerland AG
4132 MuttENZ, Switzerland

SRN AT-MF-000039061

EMDN code: W0202059085

Ref. No.	UDI-DI	Item name	Product specifications
109214	9120131862146	MS22	Low Profile Disposable Microtome Blade, angle 22°
109206	9120131862061	MS33	Low Profile Disposable Microtome Blade, angle 30°
109207	9120131862078	MS24	Low Profile Disposable Microtome Blade, angle 35°
109211	9120131862115	MS200	Low Profile Disposable Microtome Blade, angle 35°
109212	9120131862122	MS300	Low Profile Disposable Microtome Blade, angle 35°
109213	9120131862139	HP	High Profile Disposable Microtome Blade, angle 35°
109220	9120131862207	130mm Trimming Knife	Sectioning of hard tissues with fibers, blade length 130mm
109221	9120131862214	260mm Trimming Knife	Sectioning of hard tissues with fibers, blade length 260mm

Ref. No.	UDI-DI	Item name	Product specifications
110888	9120131868889	Routine LP 35	Low Profile Disposable Microtome Blade, angle 35°
110999	9120131869992	Routine HP 35	High Profile Disposable Microtome Blade, angle 35°

Intended use

Disposable, non-sterile in vitro diagnostic device for clinical use designed to be used to cut biological specimens previously inserted in paraffin. The device has been designed for professional use only in medical and veterinary laboratories and medical facilities by skilled, qualified and trained personnel. The device enables, within the workflow of a pathology laboratory, which also includes other IVDs, to make a diagnosis on a histological specimen sectioned for research and diagnostic purposes. The device is designed for use together with other devices (microtome, cryostat and paraffin) and is intended for use on the microtome and/or cryostat only and cannot be used for tissue reduction operations. The device cannot be reused; once worn out, it cannot be recycled.

Information concerning interfering substances that may affect the performance of the device

Not applicable

List of materials provided and of specific materials that are necessary but not provided

The device is supplied in a dispenser containing 50 pieces.

Information on the combined use of the device

The device is only designed for reliable, safe use together with the microtome and/or cryostat.

The device is designed to be compatible with other devices with which it comes into contact, namely paraffin and tissues subjected to cutting, but cannot be used manually as is for tissue reduction.

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Functional characteristics

- Stainless steel-carbon alloy
- Processed by way of ionic cleaning and coated with an exclusive patented PTFE film that minimizes tissue friction and compression, maintaining high cutting precision.
- The PTFE coating provides the blade with high hardness and durability.
- Manual ABS dispenser that extracts the blade, preventing any contact with the sharp parts

Instructions for use

- The device must be handled exclusively by using the recommended personal protective equipment.
- The skilled user ejects the device from the dispenser by means of the specific lever.
- Below are a few recommendations to follow for correct handling of the device in order to minimize cutting hazards.
- Carefully read the user instructions before handling the device.

1. To extract the blade from the dispenser, slide the lever towards the left until it stops.

2. With extreme caution, insert the device into the microtome blade holder, still using the provided lever, in order to facilitate insertion. If the dispenser does not eject the blade, do not force it to come out.

3. When the device has reached its expected end of service life, it can be placed at the bottom of its dispenser and then disposed of.

It is neither a self-diagnostic device nor is it intended for decentralised analysis. This device may only be used by qualified and suitably trained personnel since it is mainly intended for health professionals and operators working in specialised medical facilities. Occupational safety guidelines must be applied in accordance with the applicable regulations. The instruments used for making the diagnosis must be suitable for diagnostic use in the laboratory. A diagnosis must only be made by authorised, qualified, and trained personnel.

Conditions for sample collection, handling, and preparation

The specimen to be cut has been previously treated with chemical reagents that reduce the possibility of contamination.

Storage

It is recommended to store the product carefully in a cool, dry place protected from humidity. It is recommended to handle the blade dispenser with extreme care, avoiding impacts.

The product is stable when properly stored in unopened packaging. Do not use after the expiry date. If the product is not stored as instructed, its yield may vary, and it will have to be tested by the operator.

Stability when handling the device

The device is intended for use in specific environments where the possibility of error due to distraction cannot be eliminated. Therefore, the device must be used under the absolute control of the skilled and qualified operator, and that all instructions herein are followed in order to minimise risks that one may face, such as cutting hazards.

It is also necessary to wear the PPE required in the workplace. For successful performance of the device, do not use it beyond its expiry date.

Instructions for disposal / Safe final disposal of the device

Specific instructions on how to safely dispose of the device are to be found in the local regulations of the country where the device is used. If the product has expired and/or is unused, it must be disposed of in accordance with local waste regulations, based on the hazard classification indicated on the label and after checking for possible contamination. In some cases, analytical evaluation may be required to correctly classify the waste and determine the hazard characteristics attributed to it.

Warnings, precautions and restrictions when using the device

- Device used exclusively for the purpose specified in *Intended Use*
- Use the device only by wearing gloves and other protective PPE.
- Wear cut-resistant gloves
- Handle the blade with care, being careful not to come into contact with the sharp parts.
- Maintain a clean work environment.
- Do not reuse the device; it is **disposable**.
- When cutting infected or radioactive material, it is advisable to change the blade with every new specimen to prevent undesired contamination.
- Follow the instructions in *Instructions for use*

Device malfunction or degradation

The device itself is subject to wear and tear due to the number of cuts performed and type of tissues to be cut. By its very nature, the device may be subject to wear and tear due to use. Should the device change its appearance (presence of rust, chipped blade, etc), do not use it.

Error	Possible cause	Possible corrective action
Inadequate storage of the device	Loss of the functional features of the blade (difficulty in ejection from the dispenser due to rapid drying of the protective film) Failed or delayed diagnosis	Storage temperature information indicated on the label. Symbol, consult the IFU and specify the storage conditions of the device in the storage paragraph
Conscientious use beyond the expiry date		Indicate appropriate information on the label regarding the storage period

Warnings from the user, with the need to report any serious incidents that have occurred involving the device

In case of a serious accident, the event involving the device in question should be reported to the manufacturer or to the competent authorities, so that the evaluation of such accidents can be more easily coordinated and appropriate safety corrective actions can be taken.

Legend of symbols used in labelling and in the instructions for use


Pictogram	Pictogram name	Explanatory notes
	Manufacturer	This indicates the manufacturer of the device.
	Expiry date	This indicates the date after which the device must not be used.
	Batch number	This indicates the manufacturer 's batch required for batch identification.
	Keep away from light	This indicates that the device needs to be protected from light sources.
	Damaged packaging	This indicates that the device must not be used if the packaging is damaged, and instructions for use has to follow.
	Store in a dry place	This indicates a medical device that must be protected from moisture.
	Temperature limit (-10°C / +50°C)	This indicates the temperature limit that the medical device can be safely exposed to.
	Consult instructions for use	This indicates that the user needs to consult the instructions.
	In vitro diagnostic medical device	This indicates a device that is intended to be used as an in vitro diagnostic medical device.
	CE Marking	The CE mark certifies that the product complies with the applicable EU regulation 2017/746.
	Unique Device Identification	This indicates a carrier containing the fixed and variable information related to the UDI.

Table of revisions

Rev.	Date	Purpose of the revision
1 st issue	30/10/2023	Compliance with the new Regulation (EU) 2017/746 (In Vitro Diagnostics Medical Devices)
2 nd issue	27/12/2023	Pictogram "expiry date" replaced, CH REP inserted
3 rd issue	26/11/2025	Supplement of Routine LP 35 and Routine HP 35

You can find this document in the reserved area of the www.micros.com website.

MICROS Produktions- und Handelsges.m.b.H • Hunnenbrunn - Gewerbezone 5 • A-9300
St. Veit/Glan Telefon +43 4212 30901 • Fax +43 4212 30901-90 • E-Mail office@micros.at

Bank: Unicredit Bank Austria, **IBAN:** AT211100 0032 7004 6000, **BIC:** BKAUATWW
Firmenbuch-Nr.: FN75375x, **Gerichtsstand:** St. Veit/Glan, **UID:** ATU15049007, **EORI:** ATEOS1000002137

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E-Mail: office@micros.at

www.micros.at

CH REP Biosystems Switzerland AG
4132 Muttens, Switzerland

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