




micros
Austria








THE VIEW OF THE ESSENTIAL




soft tissue	hard tissue	cryostat	thin sectioning	ribbon sectioning	angle	pcs/box	L x W thickness	package size
-------------	-------------	----------	-----------------	-------------------	-------	---------	-----------------	--------------


LOW PROFILE - DISPOSABLE BLADES

	MS200 (Art.-Nr. 109211) LATEST PLASMA TECHNOLOGY Universal blade for daily routine in histological laboratories, ideally suited for all types of tissue. Short trimming time, very good cutting quality. Especially good for ribbon sectioning cuts, even with very thin cuts. Also suitable for industrial research applications for cutting paraffin and plastic embedded samples. All-rounder with unique hardened cutting edge and excellent durability.	● ●	● ●	●	● ●	35°	50	80 x 8 mm 0,25 mm	116 x 53 x 25 mm
	MS300 (Art.-Nr. 109212) LATEST PLASMA TECHNOLOGY A very sharp specialized blade, equally good for soft and hard tissue also recommended for cutting bone tissue. Ideal for FNA (Fine-needle aspiration) specimens, biopsies, mamma, lymph nodes and other tissue types as well as for ribbon sectioning. Latest plasma technology.	●	● ●		● ●	35°	50	80 x 8 mm 0,25 mm	116 x 53 x 25 mm
	MS24 (Art.-Nr. 109207) Robust blade with excellent sharpness and durability. Suitable for hard tissues but also for soft tissues and bone biopsies. Even for hard and difficult tissue like uterus and colon cuts. Very good results also with animal preparations. Excellent references for cryostat. Teflon coated.	●	● ●	●	●	35°	50	80 x 8 mm 0,25 mm	116 x 53 x 25 mm
	MS33 (Art.-Nr. 109206) Blade for all types of tissue especially for applications where extremely thin, high-quality sections are required. Suitable also for skin biopsy. Short trimming time and increased durability.	● ●	●		● ●	30°	50	80 x 8 mm 0,25 mm	116 x 53 x 25 mm
	MS22 (Art.-Nr. 109214) Special blade with an angle of 22° and unique plasma coating is often used for hard tissue samples. In histopathological diagnostics and industry, it impresses with its extreme sharpness, which enables precise cuts.		● ●		● ●	22°	50	80 x 8 mm 0,25 mm	116 x 53 x 25 mm

HIGH PROFILE - DISPOSABLE BLADE

	HIGH PROFILE (Art.-Nr. 109213) LATEST PLASMA TECHNOLOGY High profile blade with excellent references for cryostat and ribbon sectioning. Special hardened cutting edge, short trimming time and excellent durability.	●	●	● ●	●	35°	50	80 x 14 mm 0,32 mm	116 x 53 x 25 mm
---	---	---	---	--------	---	-----	----	-----------------------	------------------

TRIMMING KNIVES

	130 MM (Art.-Nr. 109220) 260 MM (Art.-Nr. 109221) For cutting tissue samples without causing deformation for large tissue blocks. Suitable handles see next page.					22°	50	130 x 18 mm 0,25 mm 260 x 18 mm 0,25mm	160 x 56 x 35 mm 290 x 56 x 35 mm
--	---	--	--	--	--	-----	----	---	--------------------------------------

HANDLES



Art.-Nr. 109223 for Trimming knife 260mm
Art.-Nr. 109222 for Trimming knife 130mm
Art.-Nr. 109226 für Disposable Blade 80mm



Art.-Nr. 109225
Grip for a touchless inserting of 80mm disposable blades.
Also suitable for all other brands.



📍 MICROs Produktions- und Handelsges.m.b.H.
Hunnenbrunn – Gewerbezone 5
9300 St. Veit an der Glan,
AUSTRIA

☎ Telefon: +43 4212 30901
📠 Fax: +43 4212 30901-90
✉ E-Mail: office@micros.at

🌐 www.micros.at
📘 MICROs Austria
🌐 Micros-Austria



micros
Austria